

Stainless Steel Freezer Racks

I. APPLICATIONS

1. Storage and organization of 2" and 1" standard freezer boxes.

II. FEATURES

1. Fully functional Freezer Rack is constructed of corrosion-resistant stainless steel with slide drawer
2. Includes colored labels for drawer/sample identification
3. Includes pull-out handles for easy removal off the shelf
4. Optimize Freezer and Liquid Nitrogen Storage Space
5. Ideal for storage of 2" and 1" standard size boxes
6. Fits most standard upright freezers

III. SLIDE DRAWER FREEZER RACKS DESIGN

1. The racks are available in 4 drawers, 5 drawers and 6 drawers etc.
2. Fit for 2" and 1" standard Freezer Boxes.
3. Holds 16 boxes, 20 boxes and 24 boxes etc.
4. Has three dimensions: 564*140*225(mm), 564*140*282(mm), 564*140*335(mm)

IV. PACKAGING

1. Slide Drawer Freezer Racks: 1 piece/case.

V. DESCRIPTION

Item	Description	Packaging	Sterile
99-7124	Slide Drawer Freezer Racks, Stainless steel, 6 Drawers,	1 Piece/Case	No
99-7128	Slide Drawer Freezer Racks, Stainless steel, 7 Drawers,	1 Piece/Case	No
99-7132	Slide Drawer Freezer Racks, Stainless steel, 8 Drawers,	1 Piece/Case	No
99-2118	FrameType Freezer Racks, Stainless steel, 6 layers, 3*6,	1 Piece/Case	No
99-2120	FrameType Freezer Racks, Stainless steel, 5 layers, 4*5,	1 Piece/Case	No
99-2120	FrameType Freezer Racks, Stainless steel, 7 layers, 3*7,	1 Piece/Case	No
99-2124	FrameType Freezer Racks, Stainless steel, 6 layers, 4*6,	1 Piece/Case	No
99-2128	FrameType Freezer Racks, Stainless steel, 7 layers, 4*7,	1 Piece/Case	No
99-2132	FrameType Freezer Racks, Stainless steel, 8 layers, 4*8,	1 Piece/Case	No
99-7216	Slide Drawer Freezer Racks, Stainless steel, 4 Drawers,	1 Piece/Case	No
99-7220	Slide Drawer Freezer Racks, Stainless steel, 5 Drawers,	1 Piece/Case	No
99-7224	Slide Drawer Freezer Racks, Stainless steel, 6 Drawers,	1 Piece/Case	No
99-5212	Slide Drawer Freezer Racks, For Side Access	1 Piece/Case	No
99-5215	Slide Drawer Freezer Racks, For Side Access	1 Piece/Case	No
99-5216	Slide Drawer Freezer Racks, For Side Access	1 Piece/Case	No
99-5220	Slide Drawer Freezer Racks, For Side Access	1 Piece/Case	No
99-2212	FrameType Freezer Racks, Stainless steel, 3*4	1 Piece/Case	No

99-2215	FrameType Freezer Racks, Stainless steel, 3*5	1 Piece/Case	No
99-2216	FrameType Freezer Racks, Stainless steel, 4*4	1 Piece/Case	No
99-2220	FrameType Freezer Racks, Stainless steel, 4*5	1 Piece/Case	No
99-2224	FrameType Freezer Racks, Stainless steel, 4*6	1 Piece/Case	No
99-6206	Vertical Type Freezer Rack, Stainless steel, 1*6	2 Piece/Case	No
99-6207	Vertical Type Freezer Rack, Stainless steel, 1*7	2 Piece/Case	No
99-6208	Vertical Type Freezer Rack, Stainless steel, 1*8	2 Piece/Case	No
99-6209	Vertical Type Freezer Rack, Stainless steel, 1*9	2 Piece/Case	No
99-6210	Vertical Type Freezer Rack, Stainless steel, 1*10	2 Piece/Case	No
99-6211	Vertical Type Freezer Rack, Stainless steel, 1*11	2 Piece/Case	No
99-6212	Vertical Type Freezer Rack, Stainless steel, 1*12	2 Piece/Case	No
99-6213	Vertical Type Freezer Rack, Stainless steel, 1*13	2 Piece/Case	No

VI. STERILITY

1. N/A


VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. Common Freezer Matching


Brand/Catalog#	2''boxes Plan	2''boxes Plan	2''boxes Plan	Capacity (boxes)	1''boxes Plan	1''boxes Plan	Capacity (boxes)
SANYO MDF-U74V	12*99-7220 12*99-7224	24*99-5220	12*99-2220 12*99-2224	528	24*99-7132	24*99-2132	768
SANYO MDF-U53V	8*99-7220 8*99-7224	16*99-5220	8*99-2220 8*99-2224	352	8*99-7132 8*99-7128	16*99-2132	480
Thermo 907	24*99-7220	24*99-5220	24*99-2220	480	24*99-7128	24*99-2128	672
Thermo 705	16*99-7220	16*99-5220	16*99-2220	320	16*99-7128	16*99-2128	448
Haier DW-86L628	20*99-7216	20*99-5216	20*99-2216	320	20*99-7124	20*99-2124	480
Haier DW-86L728	24*99-7220	24*99-5220	24*99-2220	480	24*99-7128	24*99-2128	672
Eppendorf U410	15*99-7216	15*99-5216	15*99-2216	240	15*99-7124	15*99-2124	360
Eppendorf U570	25*99-7216	25*99-5216	25*99-2216	400	25*99-7124	25*99-2124	600

IX. PRODUCT PARAMETER




Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H	I	J
99-7216	Racks	564	140	225	180	180	545	44	52.1	55.8	136



Dimensions=mm


Cat #	Description	A	B	C	D	E	F	G	H	I	J
99-7220	Racks	564	141	282	180	180	545	44	52.1	55.8	136



99-1220

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H	I	J
99-7224	Racks	564	140	335	180	180	545	44	52.1	55.8	136



99-1224